



LABORATORIES

Manufacturers of the Highest Quality

## PRODUCT SPECIFICATION

### Goat antibody to Human Creatine Kinase MM (CKMM), purified

**Catalog Number:** GC139      **Part Number:** 90123

**Source:** Goat Plasma

**Immunogen:** Purified CKMM from Human Skeletal Muscle

**Purification Method:** Affinity chromatography on immobilized human CKMM

#### **Form**

Clear liquid - Buffer: 10mM Sodium Phosphate, 500mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 ± 0.1

#### **Purity Specification**

Report purity by SDS-PAGE

#### **Protein / Content Specification**

≥ 5 mg/ml by A280 (E<sub>mg/ml</sub> = 1.38)

#### **Characteristics Specification**

Report IEP results vs. CKMM, CKMB, CKBB and 2x human plasma

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812  
scrippslabs.com | contact@scrippslabs.com