



LABORATORIES

Manufacturers of the Highest Quality

## PRODUCT SPECIFICATION

### Goat antibody to Human Ferritin, purified

**Catalog Number:** GF039      **Part Number:** 90125

**Source:** Goat Serum

**Immunogen:** Purified Ferritin from Human Liver

**Purification Method:** Affinity chromatography on immobilized human ferritin

#### **Form**

Clear liquid - Buffer: 50mM Tris, 500mM Sodium Chloride, 10mM EDTA, 0.1% Sodium Azide, 7.5 ± 0.1

#### **Purity Specification**

Affinity purified antibody.

Immunoabsorbed against immobilized human serum

#### **Protein / Content Specification**

Report total protein concentration

#### **Characteristics Specification**

ELISA titer (50% maximum signal);

Single precipitin arc on IEP against human ferritin;

No reaction against 2x concentrated pooled human serum

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812  
scrippslabs.com | contact@scrippslabs.com