



PRODUCT SPECIFICATION

Goat antibody to Human Cardiac Myoglobin, IgG fraction

Catalog Number: GM075 **Part Number:** 90126

Source: Goat Serum

Immunogen: Purified Human Cardiac Myoglobin

Purification Method: Ion-Exchange Chromatography

Form

Clear liquid - Buffer: 50mM Tris, 500mM Sodium Chloride, 10mM EDTA, 0.1% Sodium Azide, pH 7.5 ± 0.1.

Purity Specification

DEAE purified IgG fraction

Protein / Content Specification

Report protein concentration

Characteristics Specification

Report ELISA results, Becker titer, and cross-reactivity with hemoglobin

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.