



## PRODUCT SPECIFICATION

### Goat antibody to Recombinant Human Prolactin, purified

**Catalog Number:** GP159      **Part Number:** 90419

**Source:** Goat Serum

**Immunogen:** Purified recombinant human prolactin

**Purification Method:** Affinity chromatography

#### **Form**

Clear liquid - Buffer: 50mM Tris, 500mM Sodium Chloride, 10mM EDTA, 0.1% Sodium Azide, pH 7.5 ± 0.1.

#### **Purity Specification**

Affinity purified antibody

#### **Protein / Content Specification**

≥ 5 mg/ml by A<sub>280</sub> (E<sub>mg/ml</sub> = 1.5)

#### **Characteristics Specification**

Report ELISA titer and cross-reactivity if available.

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*