



LABORATORIES

Manufacturers of the Highest Quality

## PRODUCT SPECIFICATION

### Immunoglobulin G (IgG), Goat

**Catalog Number:** I0514      **Part Number:** 90142

**Source:** Goat Serum

#### **Form**

White to off-white lyophilized powder

#### **Purity Specification**

Report purity

#### **Protein / Content Specification**

Report approximate % protein by A<sub>280</sub> / dry weight  
Balance is buffer salts.

#### **Characteristics Specification**

This product is bulk processed and packaged by dry weight. Protein is determined by A<sub>280</sub> (E<sub>mg/ml</sub> = 1.3).

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812  
scrippslabs.com | contact@scrippslabs.com