



PRODUCT SPECIFICATION

Immunoglobulin G (IgG), Mouse

Catalog Number: I0624 **Part Number:** 90143

Source: Mouse Serum

Form

Clear to light tan liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.4 ± 0.1

Purity Specification

≥ 95% by SDS-PAGE

Protein / Content Specification

9 - 11 mg/ml (by A_{280} , $E_{mg/ml} = 1.4$)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.