



PRODUCT SPECIFICATION

Immunoglobulin G (IgG), Rabbit

Catalog Number: I0724 **Part Number:** 90144

Source: Rabbit Serum

Form

Clear liquid - Buffer: 10mM Phosphate Buffered Saline, pH 7.2 ± 0.1. Material has been sterile filtered but is not packaged as a sterile product.

Purity Specification

≥ 99% by Serum Protein Electrophoresis

Protein / Content Specification

Report protein concentration (by A_{280} , $E_{mg/ml} = 1.4$)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.