

## **PRODUCT SPECIFICATION**

### **Human Immunoglobulin G (IgG), Plasma, purified**

**Catalog Number:** I1424      **Part Number:** 90146

**Source:** Human Plasma

**Form**

Clear liquid - Buffer: 20mM Sodium Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.4 ± 0.1

**Purity Specification**

Report purity by immunoelectrophoresis (IEP)

**Protein / Content Specification**

Report protein concentration by A280 (Emg/ml = 1.35)

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 1 week) or frozen below -15°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.