



## PRODUCT SPECIFICATION

### Human Immunoglobulin A (IgA), Plasma, purified

**Catalog Number:** I2224      **Part Number:** 90147

**Source:** Human Plasma

**Form**

Clear liquid - Buffer: 20mM Tris, 100mM Sodium Chloride, 0.1% Sodium Azide, pH 8.0 ± 0.1

**Purity Specification**

Report purity by immunoelectrophoresis (IEP)

**Protein / Content Specification**

Report protein concentration by  $A_{280}$  ( $E_{mg/ml} = 1.35$ )

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 1 week) or frozen below -15°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*