

## **PRODUCT SPECIFICATION**

### **Luteinizing Hormone, alpha subunit ( $\alpha$ -LH), Bovine**

**Catalog Number:** L2114      **Part Number:** 90456

**Source:** Bovine Pituitary Glands

**Form**

Off-white to tan lyophilized powder

**Characteristics Specification**

This product is bulk processed and packaged by dry weight.

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.