

PRODUCT SPECIFICATION

Luteinizing Hormone, beta subunit (β -LH), Bovine

Catalog Number: L2314 **Part Number:** 90457

Source: Bovine Pituitary Glands

Form

Off-white to tan lyophilized powder

Characteristics Specification

This product is bulk processed and packaged by dry weight.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.