

PRODUCT SPECIFICATION

Luteinizing Hormone (LH), Monkey

Catalog Number: L3513 **Part Number:** 90376

Source: Monkey (Rhesus) Pituitary Glands

Form

Off-white to brown powder

Activity Specification

Report activity and method

Characteristics Specification

This product is bulk processed and packaged by dry weight

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.