



## **PRODUCT SPECIFICATION**

### **Myoglobin, Human Heart, $\geq 95\%$**

**Catalog Number:** M0724      **Part Number:** 90026

**Source:** Human Heart

**Form**

Clear red-brown liquid - Buffer: 20mM Tris, 1mM EDTA, 50% Glycerol, pH 8.5  $\pm$  0.1

**Purity Specification**

$\geq 95\%$  by SDS-PAGE

**Protein / Content Specification**

$\geq 2.0$  mg/ml by Lowry Protein Assay (BSA standard)

**Storage / Shipping**

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen between -10°C to -25°C. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*