



## PRODUCT SPECIFICATION

**Myoglobin, Human Heart, ≥ 99%**

**Catalog Number:** M0725      **Part Number:** 90021

**Source:** Human Heart

**Form**

Red-brown liquid. Buffer: 20mM Tris, 1mM EDTA, 50% Glycerol, pH 8.5 ± 0.1

**Purity Specification**

≥ 99% by SDS-PAGE

**Protein / Content Specification**

2.0 to 3.0 mg/ml by the Bradford method  
(BioRad Protein Assay, BSA standard)

**Storage / Shipping**

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen between -10°C to -25°C. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*