

PRODUCT SPECIFICATION

Monoclonal antibody to Human C-Reactive Protein (CRP)

Catalog Number: MC017 **Part Number:** 90274

Clone Number: CC002

Source: Murine

Immunogen: Purified human C-Reactive Protein

Purification Method: Protein A Chromatography and gel filtration

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.4 ± 0.1

Purity Specification

Protein A purified IgG fraction.
Report purity by SDS-PAGE

Protein / Content Specification

≥ 5 mg/ml by A₂₈₀ (E_{mg/ml} = 1.4)

Characteristics Specification

Mouse IgG_{2b,λ} isotype

Storage / Shipping

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.