



PRODUCT SPECIFICATION

Monoclonal antibody to Human β -2-Microglobulin (β 2M)

Catalog Number: MM017 **Part Number:** 90191

Clone Number: CM001

Source: Murine

Immunogen: Purified human β -2-Microglobulin

Purification Method: Protein A Chromatography and gel filtration

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 \pm 0.1

Purity Specification

Protein A purified IgG fraction.
Report SDS-PAGE results.

Protein / Content Specification

≥ 5 mg/ml by A_{280} ($E_{mg/ml} = 1.4$)

Characteristics Specification

Mouse IgG₁ isotype

Storage / Shipping

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.