

PRODUCT SPECIFICATION

Monoclonal antibody to Human Thyroid Stimulating Hormone (TSH)

Catalog Number: MT015 **Part Number:** 90207

Clone Number: AT001

Source: Murine

Immunogen: Purified β -TSH from human pituitary glands

Purification Method: DEAE chromatography

Form

Clear liquid - Buffer: 50mM Phosphate Buffered Saline, 0.1% Sodium Azide, pH 7.4 \pm 0.1

Purity Specification

DEAE purified IgG fraction

Protein / Content Specification

Report protein concentration by A280 (Emg/ml = 1.4)

Characteristics Specification

$K_a = 1 \times 10^{10}$

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.