



LABORATORIES

Manufacturers of the Highest Quality

## PRODUCT SPECIFICATION

### Prostatic Acid Phosphatase, Human, Liquid

**Catalog Number:** P0524      **Part Number:** 90219

**Source:** Human Seminal Fluid

#### **Form**

Clear liquid - Buffer: 10mM Tris, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 ± 0.1

#### **Purity Specification**

≥ 98% by SDS-PAGE

#### **Protein / Content Specification**

≥ 1.0 mg/ml by Lowry Protein Assay, BSA standard

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812  
scrippslabs.com | contact@scrippslabs.com