



## **PRODUCT SPECIFICATION**

**PSA/ACT Complex, Human,  $\geq 95\%$**

**Catalog Number:** P0624      **Part Number:** 90010

**Source:** PSA from human seminal fluid; ACT from human plasma

### **Form**

Clear liquid - Buffer: 10mM Sodium Acetate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 5.6  $\pm$  0.1

### **Purity Specification**

$\geq 95\%$  by SDS-PAGE

### **Protein / Content Specification**

$\geq 1.0$  mg/ml total protein by A<sub>280</sub> ( $E_{\text{mg/ml}} = 1.0$ )

Report PSA content by immunoassay (report method)

Ratio of PSA:Total Protein = 0.3 - 0.7

**Product is packaged based on the total protein value**

### **Storage / Shipping**

Product is shipped using dry ice (*not temperature monitored*). On arrival, store frozen at -20°C in a non-cycling freezer. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

### **Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*