



L A B O R A T O R I E S

Manufacturers of the Highest Quality

## PRODUCT SPECIFICATION

### Prostate Specific Antigen (PSA), Human

**Catalog Number:** P0715      **Part Number:** 90049

**Source:** Human Seminal Fluid

#### **Form**

White powder, lyophilized from 50mM ammonium bicarbonate

#### **Purity Specification**

≥ 99% by SDS-PAGE

#### **Protein / Content Specification**

> 1.0 mg/ml prior to lyophilization (report value)

by A280 (E<sub>mg/ml</sub> = 1.86)

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812

scrippslabs.com | contact@scrippslabs.com