



LABORATORIES

Manufacturers of the Highest Quality

## PRODUCT SPECIFICATION

### Prolactin, rat

**Catalog Number:** P2812      **Part Number:** 90367

**Source:** Rat Pituitary Glands

#### **Form**

Lyophilized Powder

#### **Activity Specification**

Report activity

#### **Characteristics Specification**

This product is bulk processed and packaged by dry weight

Activity is based on the pigeon crop sac weight-gain method. American J. Physiology, 105: 191 (1933).

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812  
scrippslabs.com | contact@scrippslabs.com