

PRODUCT SPECIFICATION

Plate Coating Stabilizer, no preservative

Catalog Number: P3323 **Part Number:** 90227

Form

Clear golden liquid containing a proprietary mixture of non-toxic components, 0.2µm filtered. This product does not contain preservatives and must be handled under sterile conditions.

Characteristics Specification

Plate Coating Stabilizer is a proprietary mixture of components designed to increase the stability of immobilized antigens and antibodies. This solution is supplied ready-to-use and in many cases can replace current blocking solutions. The sterile liquid is supplied without preservatives and must be handled using aseptic techniques. Alternatively, a preservative compatible with the assay may be added.

Storage / Shipping

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.