



L A B O R A T O R I E S

Manufacturers of the Highest Quality

PRODUCT SPECIFICATION

Plate Coating Stabilizer w/proprietary preservative

Catalog Number: P3323 **Part Number:** 90388

Form

Clear golden liquid containing a proprietary mixture of non-toxic components. This product contains a proprietary preservative.

Characteristics Specification

Plate Coating Stabilizer is a proprietary mixture of components designed to increase the stability of immobilized antigens and antibodies. This solution is supplied ready-to-use and in many cases can replace current blocking solutions.

Storage / Shipping

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812
scrippslabs.com | contact@scrippslabs.com