



## PRODUCT SPECIFICATION

### Thyroid Stimulating Hormone, $\alpha$ -Subunit ( $\alpha$ -TSH), Human

**Catalog Number:** T0214      **Part Number:** 90381

**Source:** Human Pituitary Glands

**Form**

White to off-white powder

**Purity Specification**

≥ 90% by SDS-PAGE

Report contaminant levels if available

**Characteristics Specification**

This material has been bulk processed and is bottled by dry weight.

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** This product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*