



PRODUCT SPECIFICATION

Thyroxine Binding Globulin (TBG), Human, $\geq 40\%$

Catalog Number: T0413 **Part Number:** 90248

Source: Human Serum/Plasma

Form

White powder, lyophilized from 20mM ammonium bicarbonate. May contain traces of buffer salts.

Purity Specification

$\geq 40\%$ (report method)

Protein / Content Specification

Report protein concentration prior to lyophilization and protein to dry weight after lyophilization.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.