

## **PRODUCT SPECIFICATION**

### **Thyroxine Binding Globulin (TBG), Human, ≥ 98%**

**Catalog Number:** T0414      **Part Number:** 90249

**Source:** Human Serum/Plasma

**Form**

White powder, lyophilized from 20mM ammonium bicarbonate. May contain traces of buffer salts.

**Purity Specification**

≥ 98% (report method)

**Protein / Content Specification**

Report protein concentration prior to lyophilization

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.