

PRODUCT SPECIFICATION

Thyroid Stimulating Hormone (TSH), canine

Catalog Number: T0914 **Part Number:** 90252

Source: Canine Pituitary Glands

Form

White to off-white lyophilized powder

Activity Specification

≥ 2.0 units/mg

Characteristics Specification

This product is bulk processed and packaged by dry weight.

Activity is determined by chick thyroidal ³²P uptake assay J. Clin. Endocrin. Metab., 31: 331 (1970)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.