



PRODUCT SPECIFICATION

Thyroid Stimulating Hormone (TSH), bovine

Catalog Number: T1614 **Part Number:** 90618

Source: Bovine Pituitary Glands

Form

White powder, lyophilized from 50mM ammonium bicarbonate

Purity Specification

Report purity by SDS-PAGE

Protein / Content Specification

Report protein concentration prior to lyophilization (Lowry Protein Assay, BSA standard)

Activity Specification

Report protein content by ELISA.

Product is packaged based on Lowry Protein.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short-term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.