

PRODUCT SPECIFICATION

Thyroid Stimulating Hormone / Luteinizing Hormone (TSH/LH), Rat

Catalog Number: T2713 **Part Number:** 90314

Source: Rat Pituitary Glands

Form

White to off-white lyophilized powder

Activity Specification

Report TSH and LH activity

Characteristics Specification

This product is bulk processed and packaged by dry weight.

TSH activity is determined by chick thyroidal ³²P uptake assay - J. Clin. Endocrin. Metab., 31: 331 (1970)

LH activity is determined by ovarian ascorbic acid depletion assay - Endocrinology, 73: 285 (1963)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.