



PRODUCT SPECIFICATION

Troponin C, Rabbit

Catalog Number: T4924 **Part Number:** 90257

Source: Rabbit Skeletal Muscle (obtained from healthy rabbits of U.S. origin)

Form

Clear colorless liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2 ± 0.1

Purity Specification

> 95% by SDS-PAGE

Protein / Content Specification

≥ 1.0 mg/ml by A_{280} ($E_{mg/ml} = 0.159$)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.