

PRODUCT SPECIFICATION

Troponin C, Human, Recombinant, ≥ 95%

Catalog Number: T5231 **Part Number:** 90691

Source: Human Troponin C, produced in *E. coli*

Form

Clear colorless liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2 ± 0.1

Purity Specification

> 95% by SDS-PAGE

Protein / Content Specification

≥ 1.0 mg/ml by A280 (Emg/ml = 0.24)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.